

Installation Column INSS 10/680



The installation column is completely made of ground stainless steel. The washing- and autopsy-table WT 10/650 can be joined by a special connection unit. The column features on the top a smooth work-surface, a seamlessly welded basin 500 x 400 x 200 mm as well as a single level tap for cold and warm water. Another single level tap with a 3 m shower-hose and a handsprinkler is mounted to the left side. At the front a covered collecting funnel, compatible to the drainage of the wash- and autopsy table WT 10/650. Service access door at the right side. All power-connections are already pre-installed inside the column.

Technical Data:

Dimensions W/D/H:	900 x 550 x 850 mm
Material:	Stainless steel (Material No. 1.4301 – ANSI 304)
Type:	INSS 10/680
Article No.:	010.010.680

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de



Installation Column INSS 10/680



Installation column INSS 10/680 with a big collection funnel



Collection funnel of the installation column with a sieve insert

TekEquipment™

Ph: 1300 368 138

www.tekequipment.com.au

info@tekequipment.com.au

TekMed NZ™

Ph: 0800 456 217

info@tekmed.co.nz

www.tekmed.co.nz

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de



KUGEL
medical